

## Autoclave



It is used for sterilization under working steam pressure upto 15psi. There are double walled units with inner and outer chamber made of heavy gauge stainless steel 304. Die pressed cover thick S.S. lid and Gasket molded silicone rubber. provided with Digital temperatures controller with indicator & Timer , Pressure gauge, steam release valve, safety valve & paddle lifting arrangement. This instrument supplied with basket of Aluminum / stainless steel. Pressure controls by spring valves .Cord &plugs to work on 220 / 230 V Ac supply.

Dimension	Inner	30 x 30h	30 x 50h	35 x 55h	45 x 75h	55 x 75h	55 x 90h
	Outer	0 x h	0 x h	0 x h	0 x h	0 x h	0 x h
Heater		2kW	2.5kw	3 kw	4 kw	5 kw	6 kw
Capaciry							
Temprature	Range	121°C					
	Accuracy	±1°C					
Pressure		15 Psi					
Guage		Calibrated Pressure Guage					
Power Supply		AC 230 V, 50/60 Hz					

### Optional :

Digital temperature controller cum indicator with Timer.  
Low water level indication with audio alarm.

**Rombai**  
INSTRUMENTS