

BOD Incubator

FEATURE (OPTIONAL FOR GMP MODEL)

ICH Guide Lines followed.

S.S. removable to prevent over shoot.

Safety thermostat to prevent over shoot.

HRC fuse for all controls.

Bacteriological incubator are non refrigerated version for temp.

Range of ambient + 5°C to 60°C



Models	Capacity CuFt	Volume LTR.	Approx Inner W X D H (cm)	No of Shelves	Temp Range	Temp Accuracy
A	3	90	40 X 40 X 57	2	5° to 60° C	± 0.5° C
B	6	175	50 X 50 X 70	2	5° to 60° C	± 0.5° C
C	10	285	57 X 57 X 88	3	5° to 60° C	± 0.5° C
D	12	340	60 X 60 X 95	3	5° to 60° C	± 0.5° C
E	16	450	70 X 70 X 90	3 / 5	5° to 60° C	± 0.5° C