

## Horizontal Laminar Air Flow

### Features:

- HEPA Filter : Mini Plate Hepa
- Pre Filter : Optional
- Ultra Violet Light
- Operating Switch
- Meno meter
- Filtration efficiency 99.99% at 0.3 micron
- Protective fiber grill
- Lowe noise level
- Durable and easy to clean
- Spare parts available at stock
- Fully made by mild steel powder coated
- Motor & blower 1440 rpm

### Specifications:

- Working Area : Class 100
- Size : 2 X 2, 3 X 2, 4 X 2, 6 X 2
- Particle Retention : 0.3 micron
- Velocity : 90 FPM  $\pm$  20 %
- Illumination : 750-800 Lumen
- Niose Level : 60 - 65 Decibles
- Power Supply : 220 Volts
- Work Surface : 304 Stainless Steel
- Acrelic Side Walls : Optional
- Florescent Lights : Optional
- Acrelic Door : Optional
- Caster Wheels : Optional
- Gas Cock : Optional

